Design summary, change log and updates.

This document is the chronological archive of the open bubble CPAP Design. It shows the working principle of the CPAP system

Iteration 1:

This is the first proposed design. It includes industrial grade as well as readily available off-the-shelf materials.


It consists of the basic fundamental parts of the system. I.e.


- Air compressor
- Flow meter
- Humidifier
- PEEP bottle
- Nasal prong



Tests, observations and changes:

- Air compressor, alternative enclosure and hardware selection.
- Breathing circuit pressure validation: relation between PEEP submersion depth (cm), Flow rate(Lpm) and pressure (manometer reading)
- Pressure deviations between bubbling & non-bubbling humidifiers.
- Water filling and submersion depth control.





Iteration 2.0

- Gas detection (NH4, CO, CO2, Propane etc.) to identify the reason of smell inside the breathing circuit.
- Effect on humidity with variable surface area inside the humidifier to control humidification.

Iteration 2.1


Breathing circuit leakage detection technique and results

Iteration 2.2


• Load testing for 3 days & observations on stability of the in-circuit humidification irrespective of ambient temperature & humidity change.

Iteration 3:

- Comparison with Pumani and experiments on rebreathing of CO2 in Open bubble CPAP.
- Humidifier water filling mechanism with ball valve.
- Load testing with the new humidifier setup, continuous run for more than 8 hours
- Comparison between looped vs unidirectional (separate inspiratory & expiratory limb) nasal prongs usage.



Iteration 4:



- Humidifier selection revised. Medical grade commonly available bottle shaped humidifiers.
- Alternate silicone tubing for the breathing circuit. (ID-8mm, OD-12 mm)
 It has advantage of compatibility with a wide range of tube connectors.
- Oxvgen blending feature addition in CPAP.

- 3-Flow-Meter method in contrast with existing pumani 2 FM design.
- Oxygen cylinder valve comparison and operating pressure levels selection.